

# Comfy™ Deviation Finger Extender

For the treatment and prevention of ulnar and radial deviations

## DEVIATE

*Pivot/swivel mechanism adjusts using hex wrench (included) to accommodate for ulnar and radial deviation*

## ADJUSTABLE

*Patented reinforced spine bends for individualized wrist and forearm needs*

## PROGRESSIVE FINGER EXTENSION

*Two graduated and removable rolls gradually increase extension of the MCP and PIP joints*

## COMFORT

*Super soft, antimicrobial, removable Headliner or Terrycloth cover is machine washable*



**High Reimbursing  
PDAC code  
L3915 or L3916**

DF-101 (h)	Deviation Finger Extender (h= headliner)
DF-101-AS (h)	Deviation Finger Extender in Adult Small (h=headliner)
PDF-101- (S,M,L)	Pediatric Deviation Finger Extender

\*\* Should be used only in connection with a care plan by a trained health care professional.



**ComfySplints**

We Make The World Comfy  
[www.comfysplints.com](http://www.comfysplints.com)

13112 Crenshaw Blvd,  
Gardena, CA 90249

T: 800-582-5332

F: 310 353-2484

[info@lenjoymedical.com](mailto:info@lenjoymedical.com)